

Clinical Science Products, Inc. offers a wide variety of FDA approved devices for urine drug testing that provide results in minutes, are easy to operate, and convenient to store. The product line that we offer combines ease of use with accuracy and rapid initial results. The **iCassette**, **iScreen**, and the **iCup** that strictly adhere to **SAMSHA** for the following urine drug testing:

Drug Name:	Abbreviation:	Cutoff Level:
Amphetamine	AMP	1000 ng/mL
Barbiturates	BAR	300 ng/mL
Benzodiazepines	BZO	300 ng/mL
Cocaine	COC	300 ng/mL
Ecstasy	MDMA	500 ng/mL
Marijuana	THC	50 ng/mL
Methadone	MTD	300 ng/mL
Methamphetamine	mAMP	1000 ng/mL
Morphine	MOR	300 ng/mL
Opiates	OPI	2000 ng/mL
Oxycodone	OXY	100 ng/mL
Phencyclidine	PCP	25 ng/mL
Propoxyphene	PPX	300 ng/mL
Tricyclic Antidepressants	TCA	1000 ng/mL



- FDA approved
- 99 % accuracy
- Laboratory accurate results
- Simple dip and read procedure
- Can test 1 – 10 drugs
- Option to photocopy results
- Fast and convenient
- 24 month shelf life

KIT SIZE
Available as 25 Test.

The **iScreen** Drug Screen Test is a cost effective urine drug testing kit – just dip and then read results. Laboratory accurate set to SAMSHA cutoff levels.



- FDA approved
- 99 % accuracy
- CLIA-waived panels available
- Laboratory accurate results
- Results in just 5 minutes
- Can test 1 – 12 drugs
- Option to photocopy results
- Simple procedure
- 24 month shelf life

KIT SIZE
Available as 25 Test.

Detect the presence of drugs in minutes.
Its drug testing efficiency you can count on.

The **iCassette** Drug Screen Test is a pipette style drug screen that is used by transferring the specimen from the collection container to the device by the use of a liquid transfer pipette or eye dropper that is supplied. The urine drug testing kits provide results in minutes.

PRECAUTION

- * For in vitro diagnostic use only.
- * Do not use after expiration date.
- * Device should not be reused.

STORAGE AND STABILITY

Store as packaged in the sealed pouch 2C to 30C for the duration of shelf life. Shelf life is 24 months and be kept away from direct sunlight, moisture and heat.

LAB CONFIRMATION SERVICES

Available upon request confirmatory testing for initially positive screening results can easily be arranged. Tests are performed in a licensed accredited laboratory using state of the art industry technology for results that are reliable and accurate.

